



Sterile Container System & Baskets



We know the value of life!



Thempson

Surgical Instruments

Surgical innovation Sytems

From Access to Repair

Address:

Veko Giz Plaza, Meydan Sokak No.3/43, Kat.13 Maslak, Istanbul,Turkey

Phone: +90 212 705 36 16, Fax: +90 212 705 36 36

info@thempson.com

www.thempson.com


ISO 13485:2003

ISO 9001 : 2008




- According to the production programme Thempson`s rights is reserved in production to improve, modify or delete the the products declared in this catalog.

- The weight information given in our catalog is approximate value.



The holding pins and plates for separation wire baskets or desire



The silicon mat, for placing the micro instruments.



Silicone handle has available in six different colors and detachable by user.




The drape retainer.





Spiral holder




Wire bottom, balanced placement and steam circulation safety.




The fixation clamps for the cylindrical instruments in different sizes.




1. Before sterilization the valves are in relax position.




2. During the sterilization by effect of vacuum the valves are opened completely and the dry air goes out of the container.



3. At the end the sterilization cycle the valves are opened in the reverse way by effect of vacuum and the sterilization cycle is completed.



4. At the end of the sterilization cycle the sterilizer equates the pressure level and the valves of the Barrier Model Container are closed.



5. After sterilization the valves are in relax position.




Sterile Container System

•
:


1/1 Full Size

Classic Model

Specifications Sizes	Lid (Perforated)	Lid (Perforated)	Safety Lid	Safety Lid
	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)	Box (Perforated)
580 X 280 X 100 mm 	10-010 10-011 10-012 10-013 10-014 10-015	10-040 10-041 10-042 10-043 10-044 10-045	10-150 10-151 10-152 10-153 10-154 10-155	10-180 10-181 10-182 10-183 10-184 10-185
580 X 280 X 135 mm 	10-016 10-017 10-018 10-019 10-020 10-021	10-046 10-047 10-048 10-049 10-050 10-051	10-156 10-157 10-158 10-159 10-160 10-161	10-186 10-187 10-188 10-189 10-190 10-191
580 X 280 X 150 mm 	10-022 10-023 10-024 10-025 10-026 10-027	10-052 10-053 10-054 10-055 10-056 10-057	10-162 10-163 10-164 10-165 10-166 10-167	10-192 10-193 10-194 10-195 10-196 10-197
580 X 280 X 200 mm 	10-028 10-029 10-030 10-031 10-032 10-033	10-058 10-059 10-060 10-061 10-062 10-063	10-168 10-169 10-170 10-171 10-172 10-173	10-198 10-199 10-200 10-201 10-202 10-203
580 X 280 X 260 mm 	10-034 10-035 10-036 10-037 10-038 10-039	10-064 10-065 10-066 10-067 10-068 10-069	10-174 10-175 10-176 10-177 10-178 10-179	10-204 10-205 10-206 10-207 10-208 10-209
Ps: All container measurements cover the outer size of the lids.				



CODE	OUTSIDE	INSIDE	WEIGHT
10-100	595 X 275 X 95 mm	543 X 260 X 75 mm	2.15 kg
10-101	595 X 275 X 116 mm	543 X 260 X 96 mm	2.30 kg
10-102	595 X 275 X 135 mm	543 X 260 X 115 mm	2.40 kg
10-103	595 X 275 X 185 mm	543 X 260 X 165 mm	2.75 kg
10-104	595 X 275 X 243 mm	543 X 260 X 223 mm	3.25 kg



CODE	OUTSIDE	INSIDE	WEIGHT
10-105	595 X 275 X 95 mm	543 X 260 X 90 mm	1.70 kg
10-106	595 X 275 X 116 mm	543 X 260 X 111 mm	2.05 kg
10-107	595 X 275 X 135 mm	543 X 260 X 130 mm	2.15 kg
10-108	595 X 275 X 185 mm	543 X 260 X 180 mm	2.60 kg
10-109	595 X 275 X 243 mm	543 X 260 X 238 mm	3.10 kg



Perforated Lids

Weight: 1.35 kg

10-120	
10-121	
10-122	
10-123	
10-124	
10-125	

Safety Lids

Weight: 1.10 kg

10-230	
10-231	
10-232	
10-233	
10-234	
10-235	

Removal Lids (Non Perforated)

Weight: 1.10 kg

10-130	
10-131	
10-132	
10-133	
10-134	
10-135	

























































































Sterile Container System

•


3/4 Quarter Size

Classic Model

Specifications Sizes	Lid (Perforated)	Lid (Perforated)	Safety Lid	Safety Lid
	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)	Box (Perforated)
465 X 280 X 100 mm 	10-240  10-241  10-242  10-243  10-244  10-245 	10-260  10-261  10-262  10-263  10-264  10-265 	10-350  10-351  10-352  10-353  10-354  10-355 	10-370  10-371  10-372  10-373  10-374  10-375 
465 X 280 X 135 mm 	10-246  10-247  10-248  10-249  10-250  10-251 	10-266  10-267  10-268  10-269  10-270  10-271 	10-356  10-357  10-358  10-359  10-360  10-361 	10-376  10-377  10-378  10-379  10-380  10-381 
465 X 280 X 150 mm 	10-252  10-253  10-254  10-255  10-256  10-257 	10-272  10-273  10-274  10-275  10-276  10-277 	10-362  10-363  10-364  10-365  10-366  10-367 	10-382  10-383  10-384  10-385  10-386  10-387 
Ps: All container measurements cover the outer size of the lids.				
				




CODE	OUTSIDE	INSIDE	WEIGHT
10-300	475 X 275 X 95 mm	422 X 260 X 90 mm	1.65 kg
10-301	475 X 275 X 116 mm	422 X 260 X 96 mm	1.80 kg
10-302	475 X 275 X 137 mm	422 X 260 X 117 mm	1.95 kg



CODE	OUTSIDE	INSIDE	WEIGHT
10-305	475 X 275 X 95 mm	422 X 260 X 90 mm	1.50 kg
10-306	475 X 275 X 116 mm	422 X 260 X 111 mm	1.65 kg
10-307	475 X 275 X 137 mm	422 X 260 X 132 mm	1.75 kg




Perforated Lids



Weight: 0.95 kg

10-340	
10-341	
10-342	
10-343	
10-344	
10-345	


Safety Lids



Weight: 0.90 kg

10-410	
10-411	
10-412	
10-413	
10-414	
10-415	

Removal Lids (Non Perforated)



Weight: 0.85 kg

10-320	
10-321	
10-322	
10-323	
10-324	
10-325	






Sterile Container System


1/2 Half Size

Classic Model

Specifications Sizes	Lid (Perforated)	Lid (Perforated)	Safety Lid	Safety Lid
	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)	Box (Perforated)
285 X 280 X 100 mm 	10-500 10-501 10-502 10-503 10-504 10-505	10-530 10-531 10-532 10-533 10-534 10-535	10-640 10-641 10-642 10-643 10-644 10-645	10-670 10-671 10-672 10-673 10-674 10-675
285 X 280 X 135 mm 	10-506 10-507 10-508 10-509 10-510 10-511	10-536 10-537 10-538 10-539 10-540 10-541	10-646 10-647 10-648 10-649 10-650 10-651	10-676 10-677 10-678 10-679 10-680 10-681
285 X 280 X 150 mm 	10-512 10-513 10-514 10-515 10-516 10-517	10-542 10-543 10-544 10-545 10-546 10-547	10-652 10-653 10-654 10-655 10-656 10-657	10-682 10-683 10-684 10-685 10-686 10-687
285 X 280 X 200 mm 	10-518 10-519 10-520 10-521 10-522 10-523	10-548 10-549 10-550 10-551 10-552 10-553	10-658 10-659 10-660 10-661 10-662 10-663	10-688 10-689 10-690 10-691 10-692 10-693
285 X 280 X 260 mm 	10-524 10-525 10-526 10-527 10-528 10-529	10-554 10-555 10-556 10-557 10-558 10-559	10-664 10-665 10-666 10-667 10-668 10-669	10-694 10-695 10-696 10-697 10-698 10-699
Ps: All container measurements cover the outer size of the lids.				




CODE	OUTSIDE	INSIDE	WEIGHT
10-600	306 X 272 X 95 mm	253 X 258 X 75 mm	1.25 kg
10-601	306 X 272 X 114 mm	253 X 258 X 94 mm	1.40 kg
10-602	306 X 272 X 135 mm	253 X 258 X 115 mm	1.55 kg
10-603	306 X 272 X 184 mm	253 X 258 X 164 mm	1.80 kg
10-604	306 X 272 X 241 mm	253 X 258 X 221 mm	2.10 kg



CODE	OUTSIDE	INSIDE	WEIGHT
10-605	306 X 272 X 95 mm	253 X 258 X 90 mm	1.15 kg
10-606	306 X 272 X 114 mm	253 X 258 X 109 mm	1.25 kg
10-607	306 X 272 X 135 mm	253 X 258 X 130 mm	1.35 kg
10-608	306 X 272 X 184 mm	253 X 258 X 179 mm	1.55 kg
10-609	306 X 272 X 241 mm	253 X 258 X 236 mm	2.00 kg




Perforated Lids



Weight: 0.70 kg

10-630	Grey
10-631	Yellow
10-632	Green
10-633	Blue
10-634	Red
10-635	Black


Safety Lids



Weight: 0.55 kg

10-700	Grey
10-701	Yellow
10-702	Green
10-703	Blue
10-704	Red
10-705	Black

Removal Lids (Non Perforated)



Weight: 0.60 kg

10-620	Grey
10-621	Yellow
10-622	Green
10-623	Blue
10-624	Red
10-625	Black
































Sterile Container System


XL Size

Classic Model

Specifications Sizes	Lid (Perforated)  Box (Non-Perforated)	Lid (Perforated)  Box (Perforated)
	590 X 370 X 120 mm 	10-136  10-137  10-138  10-139  10-140  10-141 
590 X 370 X 180 mm 	10-142  10-143  10-144  10-145  10-146  10-147 	10-216  10-217  10-218  10-219  10-220  10-221 
Ps: All container measurements cover the outer size of the lids.		
		



CODE	OUTSIDE	INSIDE	WEIGHT
10-400	595 X 365 X 105 mm	530 X 350 X 85 mm	2.50 kg
10-401	595 X 365 X 167 mm	530 X 350 X 145 mm	3.25 kg



CODE	OUTSIDE	INSIDE	WEIGHT
10-402	595 X 365 X 105 mm	543 X 350 X 100 mm	2.20 kg
10-403	595 X 365 X 167 mm	530 X 350 X 162 mm	3.00 kg



Perforated Lids

Weight: 1.65 kg

10-326	
10-327	
10-328	
10-329	
10-330	
10-331	






Sterile Container System


1/1 Full Size

Bio-Barrier Model

Specifications Sizes	Lid (with barrier)	Lid (with barrier)	Safety Lid Lid (with barrier)	Safety Lid Lid (with barrier)
	Box (Non-Perforated)	Box (with barrier)	Box (Non-Perforated)	Box (with barrier)
580 X 280 X h mm 	11-010 11-011 11-012 11-013 11-014 11-015 h= 115 mm	11-040 11-041 11-042 11-043 11-044 11-045 h= 135 mm	11-150 11-151 11-152 11-153 11-154 11-155 h= 128 mm	11-180 11-181 11-182 11-183 11-184 11-185 h= 145 mm
580 X 280 X h mm 	11-016 11-017 11-018 11-019 11-020 11-021 h= 140 mm	11-046 11-047 11-048 11-049 11-050 11-051 h= 163 mm	11-156 11-157 11-158 11-159 11-160 11-161 h= 153 mm	11-186 11-187 11-188 11-189 11-190 11-191 h= 170 mm
580 X 280 X h mm 	11-022 11-023 11-024 11-025 11-026 11-027 h= 155 mm	11-052 11-053 11-054 11-055 11-056 11-057 h= 175 mm	11-162 11-163 11-164 11-165 11-166 11-167 h= 168 mm	11-192 11-193 11-194 11-195 11-196 11-197 h= 185 mm
580 X 280 X h mm 	11-028 11-029 11-030 11-031 11-032 11-033 h= 205 mm	11-058 11-059 11-060 11-061 11-062 11-063 h= 235 mm	11-168 11-169 11-170 11-171 11-172 11-173 h= 218 mm	11-198 11-199 11-200 11-201 11-202 11-203 h= 245 mm
580 X 280 X h mm 	11-034 11-035 11-036 11-037 11-038 11-039 h= 265 mm	11-064 11-065 11-066 11-067 11-068 11-069 h= 295 mm	11-174 11-175 11-176 11-177 11-178 11-179 h= 278 mm	11-204 11-205 11-206 11-207 11-208 11-209 h= 305 mm
Ps: All container measurements cover the outer size of the lids.				



CODE	OUTSIDE	INSIDE	WEIGHT
11-100	595 X 275 X 115 mm	543 X 260 X 90 mm	2.70 kg
11-101	595 X 275 X 136 mm	543 X 260 X 111 mm	3.20 kg
11-102	595 X 275 X 155 mm	543 X 260 X 130 mm	3.40 kg
11-103	595 X 275 X 205 mm	543 X 260 X 180 mm	4.10 kg
11-104	595 X 275 X 263 mm	543 X 260 X 238 mm	4.90 kg



CODE	OUTSIDE	INSIDE	WEIGHT
11-105	595 X 275 X 95 mm	543 X 260 X 90 mm	1.70 kg
11-106	595 X 275 X 116 mm	543 X 260 X 111 mm	2.05 kg
11-107	595 X 275 X 135 mm	543 X 260 X 130 mm	2.15 kg
11-108	595 X 275 X 185 mm	543 X 260 X 180 mm	2.60 kg
11-109	595 X 275 X 243 mm	543 X 260 X 238 mm	3.10 kg



Lid Barrier system

Weight: 1.85 kg

11-120	
11-121	
11-122	
11-123	
11-124	
11-125	

Safety Lids

Weight: 1.21 kg

11-400	
11-401	
11-402	
11-403	
11-404	
11-405	



- GREY
- YELLOW
- GREEN
- BLUE
- RED
- BLACK

